

# Anti-CD47 antibody

Catalog Number: 141087

## Product name

Anti-CD47 antibody

## Specificity

Mouse CD47

## Antibody description

Rabbit Polyclonal to Mouse CD47

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Mouse CD47, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P


## Dilutions

IHC-P: 0.05-1µg/mL

## Validations



Mouse CD47 Immunohistochemistry(IHC) 15638



Mouse CD47 Immunohistochemistry(IHC) 15639