

Anti-MYB antibody

Catalog Number: 141135

Product name

Anti-MYB antibody

Specificity

Human MYB / c-Myb

Antibody description

Rabbit Polyclonal to Human MYB

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human MYB / c-Myb, and purified by antigen affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P


Dilutions

IHC-P: 0.1-2 μ g/mL

Validations



MYB / c-Myb Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry



MYB / c-Myb Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry