

Anti-GALNT7 antibody



Catalog Number: 141191

Product name

Anti-GALNT7 antibody

Immunogen

[Human GALNT7/GalNAc-T7 \(His Tag\) recombinant protein](#)

Specificity

Human GALNT7 / GalNAc-T7

Antibody description

Rabbit Polyclonal to Human GALNT7

Preparation

Produced in rabbits immunized with purified, recombinant Human GALNT7 / GalNAc-T7 . GALNT7 / GalNAc-T7 specific IgG was purified by Human GALNT7 / GalNAc-T7 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP


Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human GALNT7 / GalNAc-T7.

IP: 4-8 µl/mg of lysate

Validations



Items	A	B	C
Sample (whole cell lysate)	A431	Hela	Jurkat
Sample quantity	0.5 mg		
IP antibody quantity	4 µL		
Protein G agarose	15 µl of 50% Protein G Agarose		
Gel	15% SDS-PAGE reducing gel		
Primary antibody	S3-GALNT7-41c antibody at 10 µg/ml		
Secondary antibody	DyLight 800-labeled antibody to rabbit IgG (1:1), at 1:5000 dilution		

GALNT7 / GalNAc-T7 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation