

Anti-CCSAP antibody



Catalog Number: 141209

Product name

Anti-CCSAP antibody

Specificity

Human CCSAP

Antibody description

Rabbit Polyclonal to Human CCSAP

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human CCSAP, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



Items	Lanes	A
Sample (whole cell lysate)		ACHN
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		5-20 µg/ml
Secondary Antibody		DyLight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
Explanation		Developed using Odyssey imaging system. Predicted band size : 30 kDa Observed band size : 35 kDa

Human CCSAP Western blot (WB) 14932