

# Anti-FOXP2 antibody



Catalog Number: 141279

## Product name

Anti-FOXP2 antibody

## Specificity

Human FOXP2

## Antibody description

Rabbit Polyclonal to Human FOXP2

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human FOXP2, and purified by antigen affinity chromatography.

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF, IF

## Dilutions


WB: 5-20  $\mu$ g/ml

ICC/IF: 1-5  $\mu$ g/mL

## Validations



Human FOXP2 Immunofluorescence(IF) 15152



Lanes		A	B
Items			
Sample (whole cell lysate)		HeLa	HEK293
Sample Volume ( $\mu$ g/lane)		30	30
Gel	7.5% SDS-PAGE reducing gel		
Recommended Concentration	5-20 $\mu$ g/ml		
Secondary Antibody	Dylight 800 labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.		
Developed using Odyssey imaging system.			
Explanation	Predicted band size : 80 kDa Observed band size : 70 kDa		

Human FOXP2 Western blot (WB) 15151