

Anti-CD44 antibody



Catalog Number: 141538

Product name

Anti-CD44 antibody

Specificity

Mouse CD44

Antibody description

Rabbit Polyclonal to Mouse CD44

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Mouse CD44, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.1-2 µg/mL

Validations



Mouse CD44 Immunohistochemistry(IHC) 15244