

# Anti-GORASP2 antibody



Catalog Number: 141566

## Product name

Anti-GORASP2 antibody

## Specificity

Human GORASP2

## Antibody description

Rabbit Polyclonal to Human GORASP2

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human GORASP2, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P


## Dilutions

IHC-P: 0.1-2 µg/mL

## Validations



Human GORASP2 Immunohistochemistry(IHC)



15237

Human GORASP2 Immunohistochemistry(IHC)

15238