

Anti-IFNA2 antibody



Catalog Number: 141706

Product name

Anti-IFNA2 antibody

Immunogen

[Mouse IFNA2 recombinant protein](#)

Specificity

Mouse IFNA2

Antibody description

Rabbit Polyclonal to Mouse IFNA2

Preparation

Produced in rabbits immunized with purified, recombinant Mouse IFNA2. IFNA2 specific IgG was purified by Mouse IFNA2 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.1-2 µg/mL

Validations



IFNA2 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry