

# Anti-HA antibody



Catalog Number: 142209

## Product name

Anti-HA antibody

## Immunogen

[H5N1 HA \(His Tag\) recombinant protein](#)

## Specificity

Influenza A H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA

## Antibody description

Rabbit Polyclonal to H5N1 HA

## Preparation

Produced in rabbits immunized with purified, recombinant Influenza A H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA. Influenza A H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA specific IgG was purified by Influenza A H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ELISA, IHC-P, FCM, ICC/IF, IF, IP


## Dilutions

WB: 0.5-5 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Influenza A H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA.

## Validations



H5N1 HA Western blot (WB) 14708