

## General Information

### Product name:

Anti-CXCL8 antibody

### Immunogen:

[Rhesus IL-8/CXCL8 recombinant protein](#)

### Applications:

ELISA

### Specificity:

Rhesus IL-8 / Interleukin-8 / CXCL8

### Antibody description:

Rabbit Monoclonal to Rhesus CXCL8

## Properties

### Preparation:

This antibody was obtained from a rabbit immunized with purified, recombinant Rhesus IL-8 / Interleukin-8 / CXCL8 .

### Formulation:

0.2 µm filtered solution in PBS

### Storage:

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

### Clonality:

Monoclonal

**Ig Type:**Rabbit IgG

## Applications

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Rhesus IL-8 / Interleukin-8 / CXCL8.

## Validations