

# Anti-CD3D antibody



Catalog Number: 142655

## Product name

Anti-CD3D antibody

## Specificity

Mouse CD3d / CD3 delta

## Antibody description

Rabbit Polyclonal to Mouse CD3D

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Mouse CD3d / CD3 delta, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1-5 µg/ml

## Validations



Items	Lanes	A	B
Sample (whole cell lysate)		HepG2	Molt4
Sample Volume (µg/lane)		30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	1.5 µg/ml		
Secondary Antibody	Dylight 800 labeled Antibody To Rabbit IgG (1:1), at 1:5000 dilution.		
Explanation	Developed using Odyssey imaging system. Predicted band size : 19 kDa Observed band size : 19 kDa Additional bands at : 34 kDa (We are unsure as to the identity of these extra bands.)		

Mouse CD3d/CD3 delta Western blot (WB) 14661