

Anti-GSC antibody



Catalog Number: 142744

Product name

Anti-GSC antibody

IHC-P: 0.1-2 µg/mL

ICC/IF: 0.5-3 µg/mL


Specificity

Human GSC

Validations

Antibody description

Rabbit Polyclonal to Human GSC




Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human GSC, and purified by antigen affinity chromatography.

GSC Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Formulation

0.2 µm filtered solution in PBS



Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.


Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

GSC Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Clonality

Polyclonal



Ig Type

Rabbit IgG

Applications

IHC-P, ICC/IF, IF

GSC Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence

Dilutions