

# Anti-FSTL1 antibody



Catalog Number: 142750

## Product name

Anti-FSTL1 antibody

## Immunogen

[Mouse FSTL1 \(His Tag\) recombinant protein](#)

## Specificity

Mouse FSTL1

## Antibody description

Rabbit Polyclonal to Mouse FSTL1

## Preparation

Produced in rabbits immunized with purified, recombinant Mouse FSTL1 . FSTL1 specific IgG was purified by Mouse FSTL1 affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ELISA


## Dilutions

WB: 5-20 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Mouse FSTL1.

## Validations



FSTL1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot