

Anti-GOLGA2 antibody



Catalog Number: 142786

Product name

Anti-GOLGA2 antibody

Specificity

Human GOLGA2

Antibody description

Rabbit Polyclonal to Human GOLGA2

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human GOLGA2, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 0.1-2µg/mL

Validations



Human GM130/GOLGA2

Immunohistochemistry(IHC) 15622



Human GM130/GOLGA2

Immunohistochemistry(IHC) 15623