

Anti-ITGB2 antibody



Catalog Number: 142792

Product name

Rabbit IgG

Anti-ITGB2 antibody

Applications

Specificity

IHC-P, ICC/IF, IF

Human CD18/Integrin beta 2/ITGB2

Dilutions

IHC-P: 0.1-2µg/mL

Antibody description


ICC/IF: 0.5-3µg/mL

Rabbit Polyclonal to Human ITGB2

Validations

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human CD18/Integrin beta 2/ITGB2, and purified by antigen affinity chromatography.



Formulation

0.2 µm filtered solution in PBS


Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Human CD18/Integrin beta 2/ITGB2
Immunohistochemistry(IHC) 15799



Clonality

Polyclonal

Human CD18/Integrin beta 2/ITGB2
Immunofluorescence(IF) 15800

Ig Type