

Anti-GSC antibody



Catalog Number: 142800

Product name

Anti-GSC antibody

Specificity

Human GSC

Antibody description

Rabbit Polyclonal to Human GSC

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human GSC, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IF, IP


Dilutions

WB: 5-20 µg/ml


ICC/IF: 0.5-3 µg/mL

IP: 1-4 µg/mg of lysate


Validations



GSC Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



GSC Antibody, Rabbit PAb, Antigen Affinity Purified, Immunofluorescence



GSC Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation