

Anti-TPBG antibody



Catalog Number: 142805

Product name

Anti-TPBG antibody

Specificity

Human TPBG

Antibody description

Rabbit Polyclonal to Human TPBG

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human TPBG, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 0.1-2 µg/mL

Validations



TPBG Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry