

Anti-NOXRED1 antibody



Catalog Number: 142817

Product name

Anti-NOXRED1 antibody

Specificity

Human NOXRED1

Antibody description

Rabbit Polyclonal to Human NOXRED1

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human NOXRED1, and purified by antigen affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-20 μ g/ml

Validations



Lanes	A
Items	
Sample (whole cell lysate)	HT1080
Sample Volume (μ g/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	5-20 μ g/ml
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 40 kDa Observed band size : 40 kDa

Human NOXRED1 Western blot (WB) 15875