

Anti-HEXB antibody



Catalog Number: 142848

Product name

Anti-HEXB antibody

Immunogen

[Human HEXB / Hexosaminidase B \(His Tag\) recombinant protein](#)

Specificity

Human HEXB / Hexosaminidase B

Antibody description

Rabbit Polyclonal to Human HEXB

Preparation

Produced in rabbits immunized with purified, recombinant Human HEXB / Hexosaminidase B. HEXB / Hexosaminidase B specific IgG was purified by Human HEXB / Hexosaminidase B affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/ml


This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human HEXB / Hexosaminidase B.

IHC-P: 0.1-2 µg/mL

Validations



HEXB / Hexosaminidase B Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry



HEXB / Hexosaminidase B Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry