

# Anti-IL6ST antibody



Catalog Number: 142852

## Product name

Anti-IL6ST antibody

## Immunogen

[Mouse gp130/CD130 \(His Tag\) recombinant protein](#)

## Specificity

Mouse IL6ST/gp130/CD130

## Antibody description

Rabbit Monoclonal to Mouse IL6ST

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IL6ST/gp130/CD130 .

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM: 0.5-2 µg/test

## Validations



Mouse IL6ST/gp130/CD130 Flow Cytometry (FC)  
15780