

# Anti-SRPK1 antibody



Catalog Number: 142920

## Product name

Anti-SRPK1 antibody

## Immunogen

[Human SRPK1 \(His & GST Tag\) recombinant protein](#)

## Specificity

Human SRPK1

## Antibody description

Mouse Monoclonal to Human SRPK1

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human SRPK1. The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

IHC-P


## Dilutions

IHC-P: 10-30 $\mu$ g/mL

## Validations



Human SRPK1 Immunohistochemistry(IHC) 15792



Human SRPK1 Immunohistochemistry(IHC) 15793