

# Anti-MMP2 antibody



Catalog Number: 143018

## Product name

Anti-MMP2 antibody

## Immunogen

[Mouse MMP-2 \(His Tag\) recombinant protein](#)

## Specificity

Mouse MMP2 / MMP-2 / CLG4A

## Antibody description

Rabbit Polyclonal to Mouse MMP2

## Preparation

Produced in rabbits immunized with purified, recombinant Mouse MMP2 / MMP-2 / CLG4A . MMP2 / MMP-2 / CLG4A specific IgG was purified by Mouse MMP2 / MMP-2 / CLG4A affinity chromatography.

## Formulation

0.2  $\mu$ m filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions


WB: 5-20  $\mu$ g/ml

IP: 1-4  $\mu$ l/mg of lysate

## Validations



MMP2 / MMP-2 / CLG4A Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



MMP2 / MMP-2 / CLG4A Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation