

Anti-CD70 antibody



Catalog Number: 143071

Product name

Anti-CD70 antibody

Specificity

Mouse CD70 / CD27L / TNFSF7

Antibody description

Rabbit Polyclonal to Mouse CD70

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Mouse CD70 / CD27L / TNFSF7, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



Lanes	
Items	A
Sample (Membrane Lysate)	Jurkat
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	5-20 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 21 kDa Observed band size : 25 kDa

Mouse CD70/CD27L/TNFSF7 Western blot (WB)
14923