

TERF1 Polyclonal Antibody



Catalog Number: 160021

Product name

TERF1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to TERF1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 240-419 of human TERF1 (NP_003209.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - TERF1 Polyclonal Antibody

Western blot analysis of extracts of mouse eye, using TERF1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - TERF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using TERF1 antibody at dilution of 1:200 (40x lens).