

# AKT2 Polyclonal Antibody



Catalog Number: 160060

## Product name

AKT2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to AKT2

## Preparation

Antigen: A synthetic peptide of human AKT2

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



Western blot - AKT2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using AKT2 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - AKT2 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using AKT2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.