

# JAK3 Polyclonal Antibody



Catalog Number: 160146

## Product name

Rabbit IgG

JAK3 Polyclonal Antibody

## Applications

## Specificity

WB, IF

Human, Mouse, Rat

## Dilutions

## Antibody description


WB 1:500 - 1:2000

IF 1:50 - 1:100

Polyclonal antibody to JAK3

## Validations

## Preparation



Antigen: A synthetic peptide corresponding to a sequence within amino acids 750-850 of human JAK3 (NP\_000206.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Western blot - JAK3 Polyclonal Antibody

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Western blot analysis of extracts of U-87MG cells, using JAK3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

Polyclonal

## Ig Type