

PRKACA Polyclonal Antibody



Catalog Number: 160155

Product name

PRKACA Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to PRKACA

Preparation

Antigen: A synthetic peptide of human PRKACA

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Validations



Western blot - PRKACA Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using PRKACA antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.