

# GRIN2A Polyclonal Antibody



Catalog Number: 160182

## Product name

GRIN2A Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to GRIN2A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1130-1400 of human GRIN2A (NP\_000824.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:1000


IHC 1:50 - 1:200

## Validations



Western blot - GRIN2A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GRIN2A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - GRIN2A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using GRIN2A antibody at dilution of 1:100 (40x lens).