

# TLR7 Polyclonal Antibody



Catalog Number: 160201

## Product name

TLR7 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to TLR7

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 30-230 of human TLR7 (NP\_057646.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Immunohistochemistry - TLR7 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human thyroid cancer using TLR7 Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - TLR7 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using TLR7 Antibody at dilution of



1:100 (40x lens).

Immunofluorescence - TLR7 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using TLR7 antibody . Blue: DAPI for nuclear staining.