

# CHRNE Polyclonal Antibody



Catalog Number: 160245

## Product name

CHRNE Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to CHRNE

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of human CHRNE (NP\_000071.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:1000 - 1:4000

## Validations



Western blot - CHRNE Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CHRNE antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.