

ARPC1A Polyclonal Antibody



Catalog Number: 160589

General Information

Product name:

ARPC1A Polyclonal Antibody

Immunogen:

Applications:

WB

Specificity:

Human

Antibody description:

Polyclonal antibody to ARPC1A

Properties

Preparation:

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-370 of human ARPC1A (NP_006400.2).

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:


Polyclonal

Ig Type:Rabbit IgG

Applications

WB 1:500 - 1:2000

Validations



Western blot - ARPC1A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ARPC1A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.