

Hydroxyl-Histone H2A-Y39 Monoclonal Antibody



Catalog Number: 160668

Product name

Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Monoclonal antibody to Hydroxyl-Histone H2A-Y39

Preparation

Antigen: A specific peptide of human Hydroxyl-Histone H2A-Y39

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:500 - 1:1000

Validations




Western blot - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using Hydroxyl-Histone H2A-Y39 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in



Immunohistochemistry - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded mouse colon using Hydroxyl-Histone H2A-Y39 antibody at dilution of 1:100 (40x lens).




Hydroxyl-Histone H2A-Y39 Monoclonal Antibody



Catalog Number: 160668


Immunohistochemistry - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded human tonsil using Hydroxyl-Histone H2A-Y39 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using Hydroxyl-Histone H2A-Y39 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded mouse testis using Hydroxyl-Histone H2A-Y39 antibody at dilution of 1:100 (40x lens).