

# Hydroxyl-Histone H2A-Y39 Monoclonal Antibody



Catalog Number: 160668

## Product name

Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Monoclonal antibody to Hydroxyl-Histone H2A-Y39

## Preparation

Antigen: A specific peptide of human Hydroxyl-Histone H2A-Y39

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:500 - 1:1000

## Validations



Western blot - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody


Western blot analysis of extracts of various cell lines, using Hydroxyl-Histone H2A-Y39 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in



TBST.

Immunohistochemistry - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded mouse colon using Hydroxyl-Histone H2A-Y39 antibody at dilution of 1:100 (40x lens).



# Hydroxyl-Histone H2A-Y39 Monoclonal Antibody




Catalog Number: 160668

---


Immunohistochemistry - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded human tonsil using Hydroxyl-Histone H2A-Y39 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using Hydroxyl-Histone H2A-Y39 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - Hydroxyl-Histone H2A-Y39 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded mouse testis using Hydroxyl-Histone H2A-Y39 antibody at dilution of 1:100 (40x lens).