

# PAK4 Polyclonal Antibody



Catalog Number: 160792

## Product name

PAK4 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to PAK4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 120-280 of human PAK4 (NP\_001014831.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - PAK4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PAK4 antibody at 1:5000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.