

MAP4K4 Polyclonal Antibody



Catalog Number: 160810

Product name

MAP4K4 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to MAP4K4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 480-650 of human MAP4K4 (NP_663719.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

Validations



Western blot - MAP4K4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MAP4K4 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.