

ApoE Polyclonal Antibody



Catalog Number: 161068

Product name

ApoE Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ApoE

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ApoE (NP_000032.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions

IHC 1:50 - 1:100

Validations



Immunohistochemistry - ApoE Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using ApoE antibody at dilution of 1:100 (40x lens).