

RIPK2 Polyclonal Antibody



Catalog Number: 161597

Product name

RIPK2 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to RIPK2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-540 of human RIPK2 (NP_003812.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - RIPK2 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using RIPK2 antibody . Blue: DAPI for nuclear staining.