

ANGPTL4 Polyclonal Antibody



Catalog Number: 161625

Product name

ANGPTL4 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ANGPTL4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 26-406 of human ANGPTL4 (NP_647475.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - ANGPTL4 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ANGPTL4 antibody . Blue: DAPI for nuclear staining.