

IGFBP5 Polyclonal Antibody



Catalog Number: 161835

Product name

IGFBP5 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to IGFBP5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 21-272 of human IGFBP5 (NP_000590.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - IGFBP5 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using IGFBP5 antibody .