

ACADM Polyclonal Antibody



Catalog Number: 161861

Product name

ACADM Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to ACADM

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 26-421 of human ACADM (NP_000007.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

Validations



Immunohistochemistry - ACADM Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney using ACADM Antibody at dilution of 1:200 (40x lens).