

GRIN1 Polyclonal Antibody



Catalog Number: 161880

Product name

GRIN1 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to GRIN1

Preparation

Antigen: A synthetic peptide of human GRIN1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Immunohistochemistry - GRIN1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using GRIN1 Antibody at dilution of 1:100



(40x lens).

Immunofluorescence - GRIN1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using GRIN1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.