

IFIT3 Polyclonal Antibody



Catalog Number: 161888

General Information

Product name:

IFIT3 Polyclonal Antibody

Immunogen:

Applications:

WB

Specificity:

Human, Mouse, Rat

Antibody description:

Polyclonal antibody to IFIT3

Properties

Preparation:

Antigen: A synthetic peptide of mouse IFIT3

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:


Polyclonal

Ig Type:Rabbit IgG

Applications

WB 1:500 - 1:2000

Validations



Western blot - IFIT3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using IFIT3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.