

MSLN Polyclonal Antibody



Catalog Number: 161922

Product name

MSLN Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:100

Specificity

Human, Rat

Antibody description

Polyclonal antibody to MSLN

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 37-296 of human MSLN (NP_001170826.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC, IF


Dilutions

Validations



Immunohistochemistry - MSLN Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using MSLN Antibody at dilution of 1:100



(40x lens).

Immunofluorescence - MSLN Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using MSLN antibody . Blue: DAPI for nuclear staining.