

PHC3 Polyclonal Antibody



Catalog Number: 161945

Product name

PHC3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to PHC3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human PHC3 (NP_079223.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - PHC3 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PHC3 antibody . Blue: DAPI for nuclear staining.