

FTO Polyclonal Antibody



Catalog Number: 162063

Product name

FTO Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to FTO

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 416-505 of human FTO (NP_001073901.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - FTO Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using FTO antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - FTO Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat testis using FTO antibody at dilution of 1:100 (40x



Immunohistochemistry - FTO Polyclonal Antibody

FTO Polyclonal Antibody



Catalog Number: 162063

Immunohistochemistry of paraffin-embedded

mouse testis using FTO antibody at dilution of 1:100 (40x lens).