

MUC2 Polyclonal Antibody



Catalog Number: 162202

Product name

MUC2 Polyclonal Antibody

WB 1:500 - 1:2000

IF 1:50 - 1:200

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MUC2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 4850-4950 of human MUC2 (NP_002448.4).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IF


Dilutions

Validations



Immunofluorescence - MUC2 Polyclonal Antibody

Immunofluorescence analysis of human colon carcinoma cells using MUC2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence - MUC2 Polyclonal Antibody

Immunofluorescence analysis of rat lung cells using MUC2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.