

NR1D1 Polyclonal Antibody



Catalog Number: 162255

Product name

NR1D1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to NR1D1

Preparation

Antigen: Recombinant protein of human NR1D1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB 1:500 - 1:2000

IP 1:20 - 1:50

Validations



Western blot - NR1D1 Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using NR1D1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.