

SLC22A5 Polyclonal Antibody



Catalog Number: 162257

Product name

SLC22A5 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SLC22A5

Preparation

Antigen: Recombinant protein of human SLC22A5

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Immunohistochemistry - SLC22A5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using SLC22A5 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - SLC22A5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using SLC22A5 antibody at dilution of 1:100 (40x lens).