

TXK Polyclonal Antibody



Catalog Number: 162260

Product name

TXK Polyclonal Antibody

Specificity

Rat

Antibody description

Polyclonal antibody to TXK

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human TXK (NP_003319.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:500 - 1:2000

Validations

Rat thymus



Western blot - TXK Polyclonal Antibody

Western blot analysis of extracts of rat thymus, using TXK antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 10s.