

TYROBP Polyclonal Antibody



Catalog Number: 162272

Product name

TYROBP Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to TYROBP

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human TYROBP (NP_003323.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:500 - 1:2000

Validations

Mouse pancreas



Western blot - TYROBP Polyclonal Antibody

Western blot analysis of extracts of mouse pancreas, using TYROBP antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.